

B:FORM[®]


**2 YEARS
WARRANTY**

CE

MADE IN
ITALY

THERMOFORMING MACHINES

FOLLOW US:



Powerful and silent vacuum pump of latest generation.
Infrared quartz resistance with fast response without preheating.
Energy saving and safety device.
Compatible with square and round sheets.

Silent

Fully manual but with the same performance of the more expensive models.

Comfort

Electric drive of the locking plate for more comfort.

Deluxe

Fully automatic for maximum accuracy and precision.

**BART
MEDICAL[®]**

Via Bassa 36 - 48123 Mezzano (RA) - Italia
Tel. : 0544 523326 - Fax : 0544 523326
sales@bartmedical.com - www.bartmedical.com