

**“ PRACTICALLY
INVISIBLE ”**


THE ULTIMATE SOLUTION FOR ALIGNERS

Special material, transparent and elastic, which allows you to create a dental aligner that produces a delicate and therefore more constant force.

Great elasticity that optimizes the adaptation of the tray to the dental morphology, to the attachments and to the interdental spaces obtaining excellent results. Better adhesion allows for a more stable contact point, promoting constancy of force and therefore more precise movement of the teeth.

Greater comfort and easier insertion and removal. The aligners built with GHOST ALIGNER do not interfere with aesthetics and phonetics and, moreover, do not cause any discomfort to the patient's gums and teeth.



 **Square
125x125 mm**

 **Round
Ø 120 mm**

 **Round
Ø 125 mm**

SPECIFICATIONS

Coated on both sides with a very thin protective film (20 microns) which can be removed immediately after thermoforming, in this way the sheet maintains its original transparency and shine avoiding possible opacities caused by the porous surface of the model.

Available in all standard formats: square (125 x 125 mm), round (Ø 120 mm and Ø 125 mm) and in thicknesses of 0.50 mm, 0.75 mm and 1 mm.

Pack of 20 sheets bagged together in a triple layer foil envelope to protect them from light and humidity.